

SonoScape

EG-P550N

Ultra-Slim Gastroscope

For transnasal endoscopy



Coming soon



**A distal end diameter
as slim as 5.4 mm**
*makes it easier for both Healthcare
Professionals and patients*



**Large instrument
channel of 2.2 mm**
*to achieve an adequate
suction capacity*



**Wide field of view of 140°
dual light guides**
*resulting in a brighter and better
endoscopic observation*

EG-P550N

Field of view	140°
Depth of field	3 – 100 mm
Angulation	Up 210° / Down 90°
	Right 100° / Left 100°
Distal-end outer diameter	5.4 mm
Insertion tube outer diameter	5.8 mm
Instrument channel inner diameter	2.2 mm
Working length	1100 mm

ISO 13485 **CE** ₀₁₉₇

E-EGP550N20230113



SonoScape Medical Corp.
2F, 12th Building, Shenzhen Software Park Phase II,
Keji Middle 2nd Road, Shenzhen, China
Tel: +86-755-26722890 Fax: +86-755-26722850
Email: market@sonoscape.net
www.sonoscape.com