

EXPERIENCE THE EFFICIENCY

ED-5GT Series,

Efficient Design, Efficient Visualization, Efficient Safety



Efficient Design

- 15 degree retroview approach with effective maneuverability to access Papilla sustainably.
- 1:1 torque transmission in elevator for effective cannulation.
- Smooth 4.2 mm channel for easy exchange of devices.
- Ergonomically designed control section with easy access to elevator knob.



Efficient Visualization

- HD CMOS technology provides clear endoscopy vision.
- Optical enhancement with VIST helps in locating Papilla earlier.
- Wide field of view to reduce the procedure time and effort.



Efficient Safety

Detachable distal end cap facilitates efficient cleaning and minimizes the chance of contamination making the procedure safer for the patients.



ED-5GT Series

	ED-5GT
Field of view	100°
Direction of view	Backward side viewing 15°
Depth of view	4-60mm
Distal-end diameter	13.5mm
Insertion tube diameter	11.3mm
Instrument channel diameter	4.2mm
Angulation	Up/Down 120°/90°, Left/Right 90°/100°
Working length	1250mm
Compatibility	HD-550, HD-500



Sonoscape Medical Corp.
2F, 12th Building, Shenzhen Software Park Phase II,
Keji Middle 2nd Road, Shenzhen, China
Tel: +86-755-26722890 Fax: +86-755-26722850
Email: market@sonoscape.net
www.sonoscape.com

ISO 13485 CE 0197